

ROENTOROLL RA⁺

Application

Suitable for all transport machines. Especially suitable for fast-running machines (RA-process) up to 35 seconds, for processing mammography films by 90 s-process and for use in „low rate“ process.

Characteristics

All the standard rapid-access films are developed in RA-process with a good speed yield and medium contrast. Standard films and films with modern crystal technology reach a high contrast level with rapid processing (90 s) in **Roentoroll RA⁺**. The concentrates will keep for at least 24 months at storage temperatures between 5°C and 30°C.

Pack sizes

Roentoroll RA⁺ is a three-part developer concentrate and is supplied in a pack for 2 x 20 litres of ready-to-use developer. Usable with all TETENAL chemical mixers. With a TETENAL adapter can also be made up with any other chemical on the market.

Developer temperature

The developing temperature depends on the film material used and on the time the X-ray film stays in the developer, which varies depending on the machine type.

The standard figures are:

35°C - 36°C	with rapid-access (RA) processing
32°C - 34°C	with 90 s total processing time

Replenishment


Roentoroll RA⁺ should be replenished at 40 to 50 ml per 35 x 35 cm sheet.

Mixing instructions

Never mix Part 1, 2 and 3 concentrates!

Use 14 litres of water. Parts 1, 2 and 3 are stirred in and produce 20 litres of ready-to-use replenisher.

For fresh mixes in the developer tank, add 20 ml of TETENAL starter (article no. 105445) to the replenisher solution per litre tank volume.



Certification

Tetenal operates a fully certified quality management system to DIN EN ISO9001.

This product complies with the medical devices directive 93/42/EEC concerning medical products (CE-marking).