

Offset Plate Developer

PLATE STAR 1+9

Utilisation

Developer concentrate for pre-coated positive offset plates. Suitable for machine, dip tank and manual developer. The mixture can be adapted to the type of plate by variable dilution.

In areas with a hard water supply, it may be necessary to dilute the concentrate with demineralized water.

Usage

	dilution	temperature	developing time	recommended replenishment
machine	1+9 to 1+12	22 – 24°C	20 – 30 s	mixture 1+6 with conc. ⁽¹⁾
dip tank	1+9 to 1+12	22 – 24°C	30 – 50 s	mixture 1+6 with conc. ⁽²⁾
manually	conc. to 1+12	22 – 24°C	---	---

⁽¹⁾ Replenishment rate: with undiluted concentrate (mixture 1+6): 60 – 80 ml/m²

⁽²⁾ Refill manually with undiluted concentrate (mixture 1+6)

Packaging

PLATE STAR 1+9 Entwickler

1 litres concentrate developer in PE-can, 6 x 1 litres carton

Art.-Nr. 106507